




## dAb ImmunoTech COVID-19 (SARS-CoV-2)IgG/IgM Dual Antibody Test

(Immunochromatographic Assay)

As an aid for diagnosis of SARS-CoV-2 infection. Following FDA policy for Emergency Use Authorization during EUA processing use only



### QUICK TEST PROCEDURE

For proper results the temperatures of the test cassette and specimen diluent must be between (2°C~30°C) before starting the test. The test has four simple steps as shown in Fig.1 and Fig.2 (Back Side).

- **Step-1.** Using the enclosed capillary sampler, obtain 1 drop of fingerstick whole blood specimen.
- **Step-2.** Apply 1 drop of the whole blood specimen to the well of Test cassette.
- **Step-3.** Add 2 drops of Specimen Diluent to the specimen well.
- **Step-4.** Wait for 10 minutes, and read the results by in viewing the detection window.

### Assay Workflow-Mini-Manual

As an aid for diagnosis of SARS-CoV-2 infection. Following FDA policy for Emergency Use Authorization during EUA processing use only




Fig. 1. Schematic illustration of testing procedure for COVID-19 IgG/IgM Dual Antibody Kit




Fig.2. The three type of typical results are shown throughout assay

Tel: 470-5661748  
order@dabimmunotech.com  
www.dabimmunotech.com